2016/11/14 (Mon.)

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“Laboratory Work in Physics B4” -- Handout No.6b

Your name (_____) Nov. 14th, 2016 (Mon.)

1. Mechanism of the lock-in amplifier

(1) Consider to pick up a sinusoidal signal of a certain frequency from a voltage with noise. It is generally difficult to remove noise by using ordinary filters when the frequencies of noise and signal are close.

(2) A lock-in amplifier use the _____ method to pick up the component of a certain frequency.

○ Assume that the frequency of the signal you want to pick up is f_{ref} . Also assume that the input signal is a Fourier sum of components with various frequencies. Now let's consider the component of frequency f :

$$V_f \sin(2\pi f t + \delta_f)$$

○ First the lock-in amplifier multiplies $A_{ref} \sin(2\pi f_{ref} t + \delta_{ref})$ to the input signal. Then the component of f changes to

$$V_f \sin(2\pi f t + \delta_f) \times A_{ref} \sin(2\pi f_{ref} t + \delta_{ref}) = \underline{\hspace{10cm}}$$

This new waveform can be further separated to sinusoidal waves as:

○ Then this new waveform is passed to a _____ filter, which only pass components with frequency lower than a cut-off frequency. If the cut-off frequency is low enough, the filter only passes

$$V_x^{PSD} = \underline{\hspace{10cm}}$$

○ The lock-in amplifier also multiplies $A_{ref} \cos(2\pi f_{ref} t + \delta_{ref})$ to the input signal. Then the f -component will become

$$V_f \sin(2\pi f t + \delta_f) \times A_{ref} \cos(2\pi f_{ref} t + \delta_{ref}) = \underline{\hspace{10cm}}$$

$$= \underline{\hspace{10cm}}$$



After passing this to the filter, the component that remains is

$$V_y^{\text{PSD}} = \underline{\hspace{15cm}}$$

○ Even when the input signal consists of multiple frequency components, one can pick up the component of frequency f_{ref} by using the _____ method, with the information about the phase shift.

(3) Details of the lock-in amplifier

$$V = V_x \sin(\omega t + \delta_{\text{ref}}) + V_y \cos(\omega t + \delta_{\text{ref}}) \text{ とした時の}$$

V_x and V_y correspond to sin and cos components separated by $2\pi f t + \delta_{\text{ref}}$

- ✓ A lock-in amplifier shows $V_x = V_{f=\text{ref}} \cos(\delta_f - \delta_{\text{ref}})$ and $V_y = V_{f=\text{ref}} \sin(\delta_f - \delta_{\text{ref}})$ using values such as A_{ref} , V_x^{PSD} , and V_y^{PSD} (V_x and V_y are shown in the left and the central windows respectively)
- ✓ V_x is the amplitude of the component without phase shift compared with $A_{\text{ref}} \sin(2\pi f_{\text{ref}} t + \delta_{\text{ref}})$. V_y is the amplitude of the component with 90-deg phase shift compared to $A_{\text{ref}} \sin(2\pi f_{\text{ref}} t + \delta_{\text{ref}})$. Actually shown are the root-mean-square (rms) values ($\sqrt{1/2}$ times the amplitude).
- ✓ In the lock-in amplifier, the parameter “Phase” corresponds to δ_{ref} . In other words, by changing the parameter Phase, we can change the “origin” of the time.
- ✓ By pressing “Auto Phase”, the amplifier chooses δ_{ref} so that $\delta_f - \delta_{\text{ref}} = 0$.
- ✓ Most lock-in amplifiers can output sinusoidal voltage wave with frequency f_{ref} : $V_{\text{out}} \sin(2\pi f_{\text{ref}} t)$. Notice that the sin has no phase factor δ_{ref} .
- ✓ The cutoff frequency of the low-pass filter is determined by the parameter “Time Constant τ ” i.e. $f_{\text{cutoff}} = 1/(2\pi\tau)$
 - Roughly speaking, τ corresponds to an integration time. If one sets τ longer, one can reduce noise but on the other hand longer time is required to until the value settle the correct value. For our purpose, 300 msec is appropriate.
- ✓ By changing the parameter “Harmonic”, we can pick-up n -th harmonic component.

☆ This lock-in technique is very commonly used for such as differential resistivity measurements and optical measurements (on and off the light with a “chopper” to produce ac signal).