2016/11/14 (Mon.)

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“Laboratory Work in Physics B4” -- Handout No.6

Your name (_____) Nov. 14th, 2016 (Mon.)

1. Preparation for measurement

Check measurement instruments

- Check instruments, and also connection among them

	manufacturer	model number
DC current/voltage source (for current application to the sample)		
Multimeter (for sample voltage measurement)		
Multimeter (for _____ thermometer measurement)		
Lock-in amplifier (for ac susceptibility measurement)		

(Connection → Write down in the lab notebook)

Check the measurement box

- Open the measurement box, and check how they are connected.

Sample preparation

- Check that the sample is not degraded by measuring resistance among pins. Make a table like below.

	I+	I-	V+	V-
I-				
V+				
V-				

- Fix the resistance sample to the probe. Insert the IC pin to a white connector of the probe.
- Put the susceptibility sample to the coil. Be careful so that the sample is only in one of the pick-up coils. The Cu wire should be fixed to the probe with screw to have good thermalization.
- Connect the rock-in amplifier to the probe.
- Move the switch of the measurement box to “standard resistance”. Output voltage with appropriate amplitude and frequency. Press “Auto Phase” button (→ Fix the origin to the phase origin of the electric current). Record the values of phase, V_x , and V_y .
- Move the switch to the “coil” side. Record the values of V_x , V_y .
- Insert a magnetic material (sowing needle) to the other side of the pick-up coil. Record the changes in V_x , V_y .
- Screw the cover case of the probe.
- Remove the measurement cable. Insert the probe to the glass dewar. Do not break the dewar.
- Close the flange with screws. Tighten screws equally.
- Put the measurement cable back to the probe.
- Get ready for all other connections.

Check resistivity

- Apply DC current and measure the sample voltage.
- Change the sign of the current. Measure the voltage.
- Repeat this for various values of the current. Verify that $(V_+ - V_-)/2$ follows Ohm's law (write down in the lab notebook)

[illegible]

Check the susceptometer

- Move the switch to the “standard resistance” side. Set the voltage and frequency. Press “Auto Phase”. Record the values of Phase, V_x , and V_y
- Move the switch to the “coil” side. Record V_x , and V_y .
- Repeat this process for various values of the output voltage.

(Write down this table to the notebook)

Frequency: _____

Output	V_x (std. res.)	V_y (std. res.)	Phase	V_x (coil)	V_y (coil)

Start measurement

- Set output current. Be sure to press “output”. Use current within the range where Ohm’s law is confirmed above. **Record the value of the output current (not saved in a file)!**
- Set amplitude and frequency of the lock-in amplifier. Use the values checked above. **Record the values of amplitude and frequency to the notebook (not saved in a file).**
- Run the measurement program. Start data saving. Check that the data is saved in a file.
- Pour liquid nitrogen to the glass dewar. Cover the top of the dewar with towel.

Appendix -- Japanese-English list of technical words for deeper communication!

[illegible]